

Figure 4 shown on page 11 illustrates the dual retention capability of TYPE-C silica-hydride phases. In this example the non-polar molecule glyburide is eluted with reversed-phase retention only, since retention decreases with increasing percentage acetonitrile. For the polar molecule metformin, retention increases with increasing amount of acetonitrile – typical normal-phase behaviour but with an aqueous containing mobile phase. For the example of glyburide and metformin, co-elution of the two compounds would occur at 70% acetonitrile, with a reversal of elution order above this value.

The effect of temperature in ANP is opposite to reversed-phase for many solutes.

The less hydrophobic modified phases such as Cogent Diamond Hydride™, Cogent Phenyl Hydride™, Cogent Bidentate C8™ show the best performance for separation by ANP. Bonded phases which are more hydrophobic show weaker ANP characteristics. Greater ANP separation is generally achieved using acetonitrile rather than methanol.

## Typical Application Areas

ANP is particularly useful for the analysis of polar compounds and in most cases offers a preferable alternative to polar embedded or HILIC (Hydrophilic Interaction Liquid Chromatography) phases. The technique is widely used and referenced in metabolomic profiling, Natural Products and many others.

## Advantages of TYPE-C Silica Phases with Aqueous Normal-Phase

- Retains polar and hydrophilic compounds not retained by reversed-phase
- Precision run to run is unsurpassed by leading column brands
- Enhanced LC-MS sensitivity
- Better for prep chromatography due to the high volatility of mobile phases and higher yields

## Method Development Strategy for Selection of Reversed-Phase or Aqueous Normal-Phase

**STEP 1.** After installation and conditioning of the column, it is a good idea to start with a gradient run. We suggest starting with an acidified mobile phase of water as component A and acetonitrile as component B. Acidify both components with up to 0.5% formic or acetic acid. If you are not using LC-MS, TFA (up to 0.1%) is another option.

**STEP 2.** Run about 6 column volumes of the mobile phase in Step 1 at 95% water.

**STEP 3.** Set up your instrument to run a shallow gradient from 95% water to 40% water over 20 minutes for a 75mm long column. For longer or shorter columns, modify the gradient time proportionally. *This long and shallow gradient will be beneficial for determining the optimal gradient or isocratic method to run your mixture with later. For sharper peaks and less retention, run a shorter (steeper) gradient from the same starting point to end point.*

**STEP 4.** Equilibrate the column by running 100% acetonitrile for approximately 2 minutes for the 75mm long column. Adjust run time according to your column length.

**STEP 5.** Set up your instrument to run a shallow inverse gradient using the same mobile phase composition as in Step 1 to run from 90% acetonitrile to 40% acetonitrile over 20 minutes for a 75mm length column. For longer or shorter columns, modify the gradient time proportionally. *This long and shallow gradient will be beneficial for determining the optimal gradient or isocratic method to run your mixture with later. For sharper peaks and less retention, run a shorter (steeper) gradient from the same starting point to end point.*

**STEP 6.** Evaluate both gradient runs for retention time, peak shape and elution order. Since analyte retention on these columns is compound and method specific, some compounds may not retain in Step 3 (reversed-phase) and some may not retain in Step 5 (ANP). However, one column could produce an isocratic run which retains both polar and non-polar compounds.