

## Catalog Number

49302-T07-010

**Union, Connector Set, Tee**, low pressure, 1.0 mm / 0.040 inch through hole with 1/4-28 screw threads. Carbon PEEK, round with an internal "T". Use with 1/8th inch capillary tubing. ARE-Applied Research brand.



### Additional Info

These "Union" Tee's are used to connect PEEK or stainless steel tubing with the same connectors on all ends in the low pressures zones of HPLC systems. An internal "T" flow pattern is often used to split flow into two directions or to combine two streams into one.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$943.74**

Printed On 2024/05/20