
Product Sheet



Catalog Number

9502C-S08-T-M

Caps, Snap Top, with Butyl Rubber / Tefzel Septa, fitted in smooth, soft polyethylene, Natural Transparent Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



Additional Info

Soft polyethylene caps are more pliable than other Snap Caps and are easier to place on the vial. Popular choice when using a large number of caps at one time. These caps may not seal as well as other Snap Caps and should be used with sample solvents with low volatility.

Butyl Rubber with Tefzel is a lower cost Septa but inert enough for most applications.

Click <u>HERE</u> for smaller pack quantities.

MicroSolv Technology Corporation 9158 Industrial Blvd. NE, Leland, NC 28451 tel. (732) 380-8900, fax (910) 769-9435 Email: customers@mtc-usa.com Website: www.mtc-usa.com



Printed On 2024/05/20