


Mobile Phase Reservoir and Storage Bottles

Safety-Coated:

These borosilicate glass bottles are designed for safety. They are coated in an epoxy plastic on the outside that will not break if dropped, even when full of solvent. Although 1L and 2L are the most common sizes, 5L and 10L bottles are also available. These bottles have 45mm standard screw top necks. Compatible with many mobile phase filtration devices, you can safely filter directly into these bottles and cap them for storage or transportation.

These solvent bottles are made from Type-33 borosilicate glass and have white enamel graduations and write-on patches. Due to their superior chemical stability, these bottles are useful for storage and delivery of reagents and both aqueous and/or organic solvents. Bottles with volumes of 100mL to 10L are supplied.



Catalog Numbers	Description
58802-01	HPLC mobile phase reservoir bottle. Graduated borosilicate glass. 1 liter capacity. 45mm Screw top plastic coated. A solid screw cap & pouring ring included. 1 Each.
58819-05	HPLC mobile phase reservoir bottle. Graduated borosilicate glass. 2L capacity. 45mm screw top plastic coated. A solid screw cap & pouring ring included. 1 Each.
58802-05	HPLC mobile phase reservoir bottle. Graduated borosilicate glass. 5L capacity. 45mm screw top plastic coated for safety. A solid screw cap & pouring ring included. 1 Each.
58802-10	HPLC mobile phase reservoir bottle. Graduated borosilicate glass. 10L capacity. 45mm screw top plastic coated for safety. A solid screw cap & pouring ring included. 1 Each.

Non-Coated Mobile Phase Delivery Bottles:

MicroSolv Technology provides a low cost option for solvent bottles made from Type-33 borosilicate glass that are excellent for the use of HPLC and GC Solvents. Due to their superior chemical stability, these bottles are useful for storage and delivery of your reagents and solvents that are both aqueous and/or organic in nature. These bottles are also useful as medial bottles.

Supplied with white enamel markings, a solid blue cap and a drip free sealing ring. Other color caps are available.



Catalog Numbers	Description
58807-050	HPLC mobile phase reservoir and storage bottle, 50mL. 1 Each.
58807-050-CS	HPLC mobile phase reservoir and storage bottles, 50mL. 10/case.
58807-100	HPLC mobile phase reservoir and storage bottle, 100mL. . 1 Each.
58807-100-CS	HPLC mobile phase reservoir and storage bottles, 100mL. 10/case.
58807-1L	HPLC mobile phase reservoir and storage bottle, 1 Liter. 1 Each.
58807-1L-CS	HPLC mobile phase reservoir and storage bottles, 1 Liter. 10/case.
58807-250	HPLC mobile phase reservoir and storage bottle, 250mL. 1 Each.
58807-250-CS	HPLC mobile phase reservoir and storage bottles, 250mL. 10/case.
58807-2L	HPLC mobile phase reservoir and storage bottle, 2 Liter. 1 Each.
58807-2L-CS	HPLC mobile phase reservoir and storage bottles, 2 Liter. 5/case.
58807-500	HPLC mobile phase reservoir and storage bottle, 500mL. 1 Each.
58807-500-CS	HPLC mobile phase reservoir and storage bottles, 500mL. 10/case.