

Catalog Number

49302-CX-10

Union, Connector, Cross, high pressure, 0.30 mm / 0.012 inch through hole with 10-32 screw threads. PEEK, round with an internal "X" flow pattern. ARE-Applied Research brand. Bulk Price.



Additional Info

This low dead volume, "Union" is used to connect PEEK or stainless steel tubing with the same connectors on all ends in the high pressures zones of HPLC systems. A "Cross (X)" internal flow pattern is often used to split flow into two or three directions or to combine two or three streams into one.

Use with two finger tighten connectors (nuts) such as [47301-05](#) which are sold separately.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/19