

## Catalog Number

49680-02S-005

**Column Inline Filter**, for HPLC columns, replaceable filters and All-PEEK holder. Includes 25 disposable 2 um Ti filters. ARE-Applied Research brand. 5/PK.



### Additional Info

Compatible with typical flow rate analytical columns up to 4,500 psi and 5 ml/minute. *Not recommended for use with UHPLC pressure ranges.*

- Bio-inert PEEK and Titanium.
- Economical way to protect your columns and instruments.
- Re-usable filter holder. Easy to change the filter.
- Easy to install after the pump, before the injector.

The universal connections (fits all HPLC column threads and pilot lengths) can be installed with a wrench or hand tightened. When these filters are in use, and pressure starts to rise, simply disconnect (after pressure in the column is at atmospheric), dispose of and replace with a new inline filter.

A of 25, 0.188 inch, PEEK encased titanium filters are supplied but any PEEK encased filter or

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$406.25**

Printed On 2024/05/20