

Catalog Number

9502S-34C-T

Caps, Screw Top, with Silicone Rubber / PTFE Partial Pre-Slit (Y Shape) Septa, fitted in "knurled" polypropylene, Translucent Caps. For use with 9-425 thread, 12x32mm autosampler vials. Excellent for GC and



Additional Info

Partial slit, only in the rubber leaves PTFE to form an evaporation barrier making this a good choice for volatile samples when an autosampler requires a pre-slit Septa.

Y shaped partial slit allows for better "resealing" of rubber after initial injection.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$66.30

Printed On 2024/05/20