


## Catalog Number

04025-C

**Capillary, Bare Fused Silica**, 25um ID by 1 meter. Simplus Brand. 1 EA.



## Additional Info

This untreated, fused silica capillary is clad in polyimide and designed for use in Capillary Electrophoresis. Supplied on a plastic, non static spool for storage.

Made with extreme precision, these capillaries have excellent concentricity throughout the entire length of the capillary. The outer diameter (OD) is 365um. Made in USA.

Click [HERE](#) for 10 meter length and [HERE](#) for 25 meter length.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$157.54**

Printed On 2024/05/20