

## Catalog Number

40036-25P

**HPLC Column, Amide, 4um, 4.6mm ID x 250mm Length, 100A. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Amide HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of Sugars and other Analytes with amine groups,

*This column was one of the more popular column sizes and often chosen for use with standard LC systems that have ELSD or Corona CAD detectors or are working with UV absorbing nitrogen containing compounds. Chosen when more separation between closely related compounds or longer retention times is needed .*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,277.73**

Printed On 2024/05/19