

## Catalog Number

40036-SP100

**Semi Prep HPLC Column, Amide, 4um, 100A. 10mm ID x 100mm Length. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Amide "Semi" Preparative, HPLC Column is made with 4um, 100A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for Purification or "Prep" of Sugars and other Analytes with amine groups when a sample is desired to be collected.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$3,636.53**

Printed On 2024/05/20