

## Catalog Number

40060-10P-2

**HPLC Column, Diol, 4um, 2.1mm ID x 100mm Length, 100A. Cogent  
TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Diol HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of polar compounds that are difficult to retain in Reversed Phase, including polar and non-polar peptides, amino acids and other biological compounds such as metabolites. Used in Reversed Phase, Normal Phase or Aqueous Normal Phase (ANP).

*This size column provides good resolution and good run times and is a Very Popular Column used with LCMS.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,064.53**

Printed On 2024/05/20