

Catalog Number

40200-10P-2

HPLC Column, Silica-C 2.0, 2.2um, 2.1mm ID x 100mm Length, 120A.
Cogent TYPE-C Brand. 1 EA..



Additional Info

This Cogent Silica-C 2.0™ “near UHPLC”, HPLC Column is made with 2.2um, 120A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes

This column is one of the more popular 2.0 column sizes and often chosen for LCMS when increased resolution is desired.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,295.33

Printed On 2024/05/20