

Catalog Number

40208-10P-3

HPLC Column, Bidentate C8 2.o, 2.2um, 2.1mm ID x 100mm Length,
120A. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Bidentate C8 2.o HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this smaller particle, silica hydride based stationary phase provides Increased Efficiency, Excellent Selectivity and Faster Equilibration. Considered "Near UHPLC", this column can be used on any HPLC Instrument and still deliver improved efficiency.

Chosen for analysis of hydrophobic compounds that are difficult to retain in Aqueous Mobile Phases, including non-polar compounds when C18 is too retentive.

This size column provides similar run times as a 4.6mm ID column, but uses much less mobile phase, with lower flow rates. No special equipment is needed for this column to save solvent use.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,189.81

Printed On 2024/05/20