

Catalog Number

43200-15P-2

Metabolomics Development Kit, 2.o, 2x150. 1 Cogent Diamond Hydride 2.o™ & 1 Cogent Bidentate C18 2.o™. Cogent TYPE-C Brand. 1 PK.



Additional Info

This is an economical way to purchase two HPLC Columns that are often used with LCMS in Metabolomics and many other fields.

Each Column has a 2.1mm Internal Diameter and is 150mm Long, packed with 2.2um particles and 120A pore size and are manufactured with [Cogent TYPE-C™ Silica Hydride](#).

- Click [HERE](#) for More Information about Cogent Diamond Hydride™
- Click [HERE](#) for More Information about Cogent Bidentate C18™.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,525.29

Printed On 2024/05/20