

## Catalog Number

46002-01

**Tubing Cutter, Polymeric Tubing, Guillotine Style, Red. ARE-Applied Research brand. 1 EA.**



### Additional Info

This low cost tubing cutter is intended for PTFE, ETFE, PFA, FEP and other soft polymeric tubing but can be used to cut PEEK tubing, but the cut tubing should not be used in critical HPLC connections. For HPLC Grade Cuts, we suggest using our [Clean-Cut Tubing Cutter](#).

*When working with tubing it is most important to have carefully cut tube ends. This is imperative for making a good connection.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$14.45**

Printed On 2024/05/19