


## Catalog Number

47305-07-1C

**Tubing, PEEK, 0.007 inch (0.18 mm) ID, 1/16th inch (1.6 mm) OD,**  
HPLC grade, yellow stripe, Your custom continuous length tubing, priced  
per meter. ARE-Applied Research brand.



### Additional Info

You determine the length of this HPLC grade tubing that is extruded from virgin PEEK material without any manufacturing additives. The yellow stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.71 mm
- Color code: yellow stripe
- Pressure rating: 5800 psi (399.9 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/20