


## Catalog Number

47305-20-1M

**Tubing, PEEK, 0.020 inch (0.50 mm) ID, 1/16th inch (1.6 mm) OD,**  
HPLC grade, orange stripe, 1 meter roll, ARE-Applied Research brand. 1  
EA.



### Additional Info

This HPLC grade tubing is extruded from virgin PEEK material without any manufacturing additives. The orange stripe is added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.55 mm
- Color code: orange stripe
- Pressure rating: 4500 psi (310.3 bar)

Other lengths available:

- [3 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$23.28**

Printed On 2024/05/20