


Catalog Number

47305-ST-02-C

Tubing, PEEK, 0.0025 inch (0.064 mm) ID, 1/16th inch (1.6 mm) OD,
Ultra-T grade, natural color, 25 meter roll, ARE-Applied Research brand.
1 EA.



Additional Info

This Ultra-T™ grade (*premium*) tubing is extruded from virgin PEEK material without any manufacturing additives. A natural PEEK material has no color leachables.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.768 mm
- Color code: Natural (no stripe)
- Pressure rating: 6300 psi (434.4 bar)
- Supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$770.40

Printed On 2024/05/20