


## Catalog Number

47305-ST-07-3M

**Tubing, PEEK, 0.007 inch (0.18 mm) ID, 1/16th inch (1.6 mm) OD,**  
Ultra-T grade, yellow intermittent stripes, 3 meter roll, ARE-Applied  
Research brand. 1 EA.



### Additional Info

This Ultra-T™ grade (*premium*) tubing is extruded from virgin PEEK material without any manufacturing additives. The broken yellow stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.71 mm
- Color code: yellow stripes
- Pressure rating: 5800 psi (399.9 bar)

Other lengths available:

- [1 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$92.45**

Printed On 2024/05/20