

## Catalog Number

47310-25-3M

**Tubing, PEEK, 0.010 inch (0.25 mm) ID, 1/32nd inch (0.79 mm) OD,**  
general grade, natural, 3 meter roll, ARE-Applied Research brand. 1 EA.



### Additional Info

This general grade tubing is extruded from virgin PEEK material without any manufacturing additives. No color is added leaving the tubing all-natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.27 mm
- Color code: Natural
- Pressure rating: 4500 psi (310 bar)

Other lengths available:

- [1 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$104.00**

Printed On 2024/05/20