

## Catalog Number

47450-04

**Union, Column Coupler**, universal fit 0.50 mm through hole 0.020 inch (0.50 mm) ID tubing. Color coded orange. PEEK. ARE-Applied Research brand. 1 EA.



### Additional Info

This convenient column connector is used to join two HPLC columns or a guard column and an analytical column in series.

Sleeve is color coded orange to identify the tubing internal diameter and is made from polypropylene.

Tubing is all PEEK material and is made from Ultra-T tubing for precise concentricity and better performance.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$36.60**

Printed On 2024/05/19