

Catalog Number

47452-07-010

Adapter, for use with low pressure HPLC. 3/32 inch barbed connection. Polypropylene. ARE-Applied Research brand. 10/PK.



Additional Info

This adapter has a hose barb for 3/32 inch (2.4 mm) ID tubing on one side and a female, Luer Lock port on the other side.

This adapter is great when using a syringe to connect to tubing.

- Maximum Temperature: 100°C
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for bulk pricing.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$42.26

Printed On 2024/05/20