

Catalog Number

47452-32

Adapter, for use with low pressure. Male luer to barbed connection for 1/8 inch tubing. Polypropylene. ARE-Applied Research brand. Bulk Price.



Additional Info

This adapter has a male Luer Lock connection on one side and a hose barb male connection on the other side that accommodates 1/8th inch (3.2 mm) ID tubing.

This adapter is excellent when using a syringe filter or other device that has a female Luer Lock or Luer Slip port and you want to connect with 1/8th inch (3.2 mm) OD tubing.

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for economical packs of 10.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/20