


## Catalog Number

48410-30-1C

**Tubing, FEP, 0.030 inch (0.75 mm) ID, 1/16th inch (1.6 mm) OD, low pressure, Your continuous length tubing. ARE-Applied Research brand.**



### Additional Info

You determine the length of this FEP (*Fluoroethylene Propylene*) low pressure tubing, which is chosen for use in Ion Chromatography applications, but can be used in all low-pressure, laboratory applications including before and after an HPLC Pump such as with Mobile Phase Delivery and Solvent Waste Systems. Also known as Teflon FEP.

Do not use this tubing with heat over 80°C.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Wall Thickness: 0.425 mm
- Color: Natural, Translucent
- Pressure rating: 1075 psi (74.1 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/20