

## Catalog Number

48510-30-1C

**Tubing, ETFE, 0.030 inch (0.75 mm) ID, 1/16th inch (1.6 mm) OD, low pressure, Your continuous length tubing. ARE-Applied Research brand.**



### Additional Info

You determine the length of this ETFE (*Tefzel*) tubing that is extruded from virgin ETFE without any manufacturing additives. This tubing is used widely in all low-pressure, laboratory applications including before and after an HPLC Pump such as with Mobile Phase Delivery and Solvent Waste Systems. ETFE tubing is chemically inert except for some halogenated compounds and has low permeability compared with PTFE tubing.

- Internal and Outer Dimensional Tolerance: +/- 0.10mm
- Color: Natural, Translucent
- Wall Thickness: 0.425 mm
- Pressure rating: 1625 psi (bar 112.1)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/20