

## Catalog Number

49101-75C

**Performance Qualification HPLC column.** Replacement, for use with the Chemical Solutions HSQ kit. Includes a certificate of performance. Cogent brand. 1/EA



### Additional Info

This specially selected HPLC column is tested and certified by our laboratory for specific performance in performance qualification of HPLC instruments.

Replacement for the original column supplied with the MicroSolv HSQ (High Speed Qualification) Kit, it is 75mm long and 4.6mm ID. The particle size of the column phase is 5  $\mu$ m.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$920.25**

Printed On 2024/05/19