


Catalog Number

49305-07-1C

Tubing, PEEK, 0.007 inch (0.18 mm) ID, 1/16th inch (1.6 mm) OD, general grade, solid yellow, Your custom continuous length tubing, priced per meter. ARE-Applied Research brand.



Additional Info

You determine the length of this general grade tubing that is extruded from virgin PEEK material without any manufacturing additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.71 mm
- Color code: yellow
- Pressure rating: 5800 psi (399.9 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/31