


## Catalog Number

49307-504G

**Tubing, PEEK, Coiled, 0.030 (0.75 mm) ID, 1/16th inch (1.6 mm) OD.**  
1225 mm length, HPLC grade, green stripe, with tubing connectors.  
ARE-Applied Research brand. 1 EA.



### Additional Info

These Flexible, Coiled PEEK Tubes will expand and contract which will save space in places such as behind an instrument and provide a Pulse Dampener Effect for Parabolic Flow Fluids which will result in a quieter base line and also reduce band broadening due to excess tubing lengths.

This coil is 35 mm wide when resting and 200 mm when fully extended, with a total tubing length of 1225 mm, end to end.

*Type 4 FlexiTube, color coded: green stripe for 0.030 inch (0.75 mm) ID*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$59.70**

Printed On 2024/05/19