

Catalog Number

49330-GC3

Guard Column Holder Nut, PEEK, for 3.0 mm ID Guard Columns with 8 mm active length. ARE-Applied Research brand. 1 EA.



Additional Info

This Replacement, Guard Column Holder Nut is for 3.0mm ID Guard Columns. Extruded from PEEK,

Carbon is added to PEEK to achieve black color.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$31.78

Printed On 2024/05/19