

Catalog Number

49330-GC4

Guard Column Holder Insert, PEEK, for 3.0 mm ID Guard Columns with 8 mm active length. ARE-Applied Research brand. 1 EA.



Additional Info

This Replacement, Guard Column Holder Insert is for 3.0mm ID Guard Columns. Extruded from PEEK,

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$45.15

Printed On 2024/05/19