


## Catalog Number

49410-COIL-4M

**Tubing, Stainless Steel, roll, 0.010 inch (0.25 mm) ID, 1/16 inch (1.59 mm) OD, HPLC Grade, 4-meter roll. ARE-Applied Research brand. 1 EA.**



### Additional Info

This Stainless-Steel tubing has many uses in the laboratory but is designed for HPLC. Welded and annealed, they are tested according to EN 10204/3.1B. These Coils are more economical than pre-cut tubing.

- The max pressure of this tubing is 19, 700 psi / 1,358.26 bar.
- Internal and Outer Dimension Tolerances: +/- 0.05mm.
- Wall Thickness: 0.668 mm

Click [HERE](#) for Tips on Cutting SS Tubing with a hand tool.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$82.80**

Printed On 2024/05/19