

Catalog Number

49440-30

Tubing, Stainless Steel, Pre-Cut, 0.040 inch (1.016 mm) ID, 1/16 inch (1.59 mm) OD, HPLC Grade, 30cm Long. ARE-Applied Research brand.
1 EA.



Additional Info

Properly cutting tubing is one of the most difficult tasks in a chromatography laboratory. These Pre-Cut Stainless Steel Tubes are the perfect solution.

Welded and annealed, the ends are cut using Electronic Discharge (*EDC*) Technique to achieve square and burr free zero dead volume connections. They are tested according to EN 10204/3.1B. Color Bands are not included on this size of Pre-Cut Tubing.

- *Internal and Outer Dimension Tolerances: +/- 0.05mm.*
- Wall Thickness: 0.287 mm

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$22.78

Printed On 2024/05/19