

Catalog Number

49590-F116

Ferrules, stainless steel, single cone, for 1/16 inch tubing. ARE-Applied Research brand. 10/PK.



Additional Info

These single cone ferrules are made to be used with our 316L stainless steel tubing connector male nuts when connecting 1/16th inch SS tubing to columns, fittings or detectors in high pressures zones of HPLC systems.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$37.56

Printed On 2024/05/19