

Catalog Number

58022-P04-K

Syringe Filters, 4mm, PTFE, 0.22um Pore Size. Blue Polypropylene, AQ Brand. 1000/CS.



Additional Info

Advanced Quality PTFE Syringe Filters, AQ™ brand, for Pharmaceutical Grade HPLC & Dissolution Testing. Manufactured with a unique 0.22µm PTFE or Polytetrafluoroethylene (purest form of Teflon®) in a 4mm polypropylene housing, these syringe filters are made for laboratories that require consistent, high performance results with ease of use especially during high volume work. Excellent choice when laboratories cannot afford poor chromatographic results.

These Filters are Compatible with any organic HPLC solvent and designed to have strong chemical strength (will not degrade) in the presence of strong solvents, acids and bases; not recommended for 100% water or Aqueous Buffer Solutions.

Validated for many regulatory requirements, they are designed, tested and QC Certified. They comply with US FDA CFR Title 21, Part 211.92 for minimum Fiber Release and with USP <788> test for particulate matter in injectables as well as USP for "Sterile Water for Injection" and

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,381.44

Printed On 2024/05/18