


## Catalog Number

58211-20

**Needles, Blunt Tip 90 degrees**, reusable, all metal, 20 gauge x 2.0 inches, Non Sterile for Laboratory Use Only. BASIK Brand. 12/PK



### Additional Info

These non-sterile all metal needles can be used with the Luer Lock or Luer Slip, MicroSolv Chrom Syringes in the laboratory. The needle and hub is stainless steel and can be reused after cleaning.

The 90° blunt tip is good for syphoning from vials, bottles, biological systems or any small pool of liquid. These syringes can be connected to plastic tubing to act a small, low pressure "pump".

*Not recommended for piercing autosampler vial septa.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$72.66**

Printed On 2024/05/20