

## Catalog Number

58802-05-A

**Bottle, Mobile Phase Reservoir**, Glass, Amber, 5L, Screw Top, GL45 Threads. MicroSolv Brand. 1 EA.



### Additional Info

This Amber, Borosilicate Glass Bottle has a 5L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir. The small opening, large Capacity and Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MicroSolv Filter Degassers](#) or other Degassers or disposable [Filter Degassers](#). A Solid Screw Cap without a liner is included.

Larger bottles like this one are often used for making mobile phase at the start of the day and throughout the day the bottles on the instrument can be topped off with the same mobile phase minimizing issues from mobile phase preparation. This is especially useful for overnight runs before the end of the workday.

*This Bottle is often used with the [Mobile Phase Delivery System](#) and is NOT "Safety Coated".*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$904.66**

Printed On 2024/05/19