

Catalog Number

58802-10-X

Bottle, Mobile Phase Reservoir, Glass, Clear, 10L, Screw Top, GL45 Threads. MicroSolv Brand. 1 EA.



Additional Info

This Clear, Borosilicate Glass Bottle has a 10L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir. The small opening, large Capacity and Enamel Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MicroSolv Filter Degassers](#), other Degassers or disposable [Filter Degassers](#). A Solid Screw Cap without a liner and a polypropylene Drip Ring are included.

This Bottle is often used with the [Mobile Phase Delivery System](#) and is *Not "Safety Coated"*.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$650.84

Printed On 2024/05/19