

Catalog Number

58807-100-A-CS

Bottles, HPLC Mobile Phase Storage, Glass, Amber, 100ml, Screw Top, GL45 Threads. MicroSolv Brand. 10/CS.



Additional Info

These Amber, Borosilicate Glass Bottles have a 100ml "Total Fill Capacity" and are perfect as HPLC Mobile Phase Storage Bottles. The small openings and Enamel Graduations and Marking Patches makes them a convenient way to store concentrated Mobile Phase or its Ingredients.

Blue, Solid Screw Caps without liners and polypropylene Drip Rings are included. This Bottle is *Not "Safety Coated"*.

Click [HERE](#) for smaller package quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$938.18

Printed On 2024/05/19