

Catalog Number

58807-500-CS

Bottles, HPLC Mobile Phase Storage, Glass, Clear, 500ml, Screw Top, GL45 Threads. MicroSolv Brand. 10/CS.



Additional Info

These Clear, Borosilicate Glass Bottles have a 500ml "Total Fill Capacity" and are perfect as HPLC Mobile Phase Storage Bottles. The small openings and Enamel Graduations and Marking Patches makes them a way to store concentrated Mobile Phase or its Ingredients.

Blue, Solid Screw Caps without liners and polypropylene Drip Rings are included. These Bottles are *Not "Safety Coated"*.

Click [HERE](#) for smaller package quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$227.91

Printed On 2024/05/19