


Catalog Number

58820-2L-CS

Bottles, Mobile Phase Reservoir, safety, coated glass, UV cut off of 280nm, clear, 2 L, screw Top, GL45 threads. Caps sold separately. ARE-Applied Research brand. CLEARANCE. 5/CS.



Additional Info

This Clear, Plastic Coated, Borosilicate Glass Bottle has a 2L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir.

These bottles are in excellent condition and are not "seconds", damaged, returned goods or defective in any way.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$300.94

Printed On 2024/05/19