

## Catalog Number

58930-S-A

**Syringe, Chrom. Disposable, 30ml, Luer Lock, Non Sterile. MicroSolv Brand. CLEARANCE. 50/PK.**



### Additional Info

These Luer Lock Tip, *Chrom Syringes* are chosen for many uses in chromatography as well as used with syringe filters or for any lab situation needing an inert, non-reactive, plastic syringe. The Luer Lock device provides additional safety by securing the syringe with the needle or filter to prevent the devices from popping off the syringe during times of higher pressures.

Latex free, they do not contain rubber, silicone oil, styrene or DEHP that is found in other plastic syringes. These all plastic syringes are manufactured from laboratory grade polypropylene and polyethylene. More chemically resistant than syringes with rubber-tipped plungers, they are also equipped with a positive "click", safety stop to help prevent accidental spills.

These economical, disposable syringes are not individually wrapped and are bulk packed in one bag.

Click [HERE](#) for larger package quantity or click [HERE](#) for case quantity.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$65.00**

Printed On 2024/05/20