

Catalog Number

59008-15P

HPLC Column, MicroBee C8 Octyl, 10um, 4.6mm ID x 150mm Length,
125A. Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent MicroBee C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to irregular shaped silica. This makes this column an ideal Alternate in USP Methods that use the L7 Classification and Microbondapak™ or similar C8 columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is excellent for long chain fatty acids and will have less Retention than a C18 Column.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$775.46

Printed On 2024/05/19