

## Catalog Number

59300-05

**HPLC Column, MicroBee Silica, 10um, 3.9mm ID x 300mm Length, 125A. Cogent Brand. 1 EA.**  
**MOST POPULAR.**



### Additional Info

This Very Popular Cogent MicroBee Silica™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is irregular shaped silica that is not modified in any way. This makes it an ideal Alternate in older USP Methods that use the L3 Classification for columns.

Often used for hydrophilic or polar compounds in Normal Phase HPLC.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$764.59**

Printed On 2024/06/01