

Catalog Number

63100-25P-W

HPLC Column, Silica, 10um, 4.6mm ID x 250mm Length, 100A with WCS hardware. Partisil Brand. 1 EA.



Additional Info

This Partisil Silica™ Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is not modified in any way. This Column is written in many USP Methods; Partisil Silica used in older methods with the L3 Classification.

Often used for hydrophilic or polar compounds in Normal Phase HPLC, this column is supplied with the Whatman Compression Sealing System (WCS).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,549.25

Printed On 2024/05/20