

## Catalog Number

68004-SP-PG1

**Semi Prep Guard Column Cartridge, RP C4, 10um, 100A. 10mm ID x 10mm Length. Cogent Brand. 5/PK.**



### Additional Info

These Guard Column Cartridges are packed with 10µm, Cogent RP C4™ Stationary Phase. They are double end capped (*fully*), with traditional bonding of a 4 carbon ligand (Butyl) to the highest purity silica available. They are designed to protect your RP C4 Semi Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, they have Stainless Steel Frits.

*Use these Guard Column Cartridges with a 10mm ID Prep Column Holder. Click [HERE](#) for more information.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$895.91**

Printed On 2024/05/19