

Catalog Number

68024-P50-250

Process HPLC Column, RP Phenyl, 10um, 100A. 50mm ID x 250mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent RP Phenyl™ Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 10 µm particle works well with higher flow rates and produces lower back pressures than smaller particles.

The extra large diameter and very long length of this column provides very high loading capacity which is great for purification of isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl rings, making it ideal for Process HPLC.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$21,839.09

Printed On 2024/05/20