

## Catalog Number

68124-10P-2

**UHPLC Column, RP Phenyl, 1.8 $\mu$ m, 100A. 2.1mm ID x 100mm Length.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent RP Phenyl™, is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring, to the highest purity silica available and is one of our Most Popular Columns. This is one of our Most Popular Phenyl Columns.

The small, 1.8 $\mu$ m Particle Size will provide high efficiency even at high flow rates and the Peak Capacity is greater than a 3 $\mu$ m Particle. Higher back pressures than with larger particle columns are typical and a special UHPLC instrument may be required for proper results.

Choose for analysis of isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl Rings, when you are trying to reduce the run time and keep resolution the same or if you want to improve resolution of existing methods. Ideal for new method development since these columns are lower cost to run due to reduced solvent use and wastage and increased instrument throughput

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,466.33**

Printed On 2024/05/20