

## Catalog Number

68308-25P-4

**HPLC Column, RP C8 Octyl, 3um, 4.0mm ID x 250mm Length, 100A.**

Cogent Brand. 1 EA.

*CUSTOM ORDER.*



### Additional Info

This Cogent RP C8™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (fully), with traditional bonding of an 8 (Octyl) carbon silane to the highest purity silica available, making this column ideal as an Alternate Column in USP Methods that use the L7 Classification with 3um particle size stationary phases, or for new methods.

The 3µm Particle in this column makes it more efficient than a similar 5um column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is and increased back pressure and methods that use longer columns and viscous solvents may exceed the pressure limits of your instrument. Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids.

*This column is often chosen for starting points of new methods or for Substitution for USP L7 Methods.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,036.93**

Printed On 2024/05/20