

## Catalog Number

68308-PG2

**Guard Column Cartridges, RP C8 Octyl, 3 $\mu$ m, 100A. 4.0mm ID x 10mm Length. Cogent Brand. 5/PK.**



### Additional Info

These Guard Column Cartridges are packed with 3 $\mu$ m, Cogent RP C8™ (Octyl) Stationary Phase. They are double end capped (*fully*), with traditional bonding of an 8 carbon silane to the highest purity silica available.

Designed to protect your 3.0, 4.0 or 4.6mm ID, RP C8, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

*These Guard Column Cartridges must be used with a Guard Column Holder, click [HERE](#) for more information about the Holder.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$447.95**

Printed On 2024/05/20