

## Catalog Number

68322-PG2

**Guard Column Cartridges, RP PFP, 3 $\mu$ m, 100A. 4.0mm ID x 10mm Length. Cogent Brand. 5/PK.**



### Additional Info

These Guard Column Cartridges are packed with 3 $\mu$ m, Cogent RP PFP™ Pentafluorophenyl Stationary Phase. It is double end capped (*fully*), and bonded with a four carbon link to a benzene ring with 5 Fluorines on the highest purity 100A silica available.

Designed to protect your 3.0, 4.0 or 4.6mm ID, RP PFP™, Analytical Columns from particulate and contaminants that will adsorb to it. They are packed in Stainless Steel Tubes and have Stainless Steel Frits that are pressed into PEEK Frit Caps.

*These Guard Column Cartridges must be used with a Guard Column Holder, click [HERE](#) for more information about the Holder.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$447.95**

Printed On 2024/05/19